

# LVL Starter Packs Available in 3 sizes

PACKS TO SUIT EVERY BIOBANK AND SAMPLE STORAGE NEEDS

## Single Tube Reader

- › The Single Tube Reader scans the 2D codes on SAFE® tubes and other cryo tubes
- › Scans SBS and Cryo Tubes in just 5 seconds
- › Easy export in MS Excel, XML or text and scanned images can also be saved
- › Plug & play



## SBS Format Rack Reader

- › The new standard for fast and reliable SBS 2D Tube barcode decoding
- › Easy export to MS Excel, XML, text, JSON, and jpg
- › Small footprint enables a simple integration in the automatic pipetting systems
- › Available with USB Decoding software with lifelong, free upgrades
- › Plug & play



## Single tube picker and de/capper

- › Easy picking, capping and de-capping of single tubes out of a rack
- › The caps hold fast to the device and are able to be ejected with a push of the white plunger



## Semi Automated Handheld Capper/Decapper

- › For opening and closing of tubes with screw caps simultaneously
- › Available in 8 channel
- › Mobile, battery operated 8 channel for table use
- › Power-saving mode for longer operating time in wireless mode



## SAFE® 2D barcoded racked Tubes






- › Excellent sealing (IATA certified and Co<sub>2</sub> tested)
- › Automation friendly or easy manual handling
- › High performance fibre laser coding for unmatched contrast and readability
- › 2D bottom code, 1D & HR side code
- › SAFE® low leachables polypropylene



Contact the Pathtech team

Phone 1800 069 161 | Email [info@pathtech.com.au](mailto:info@pathtech.com.au) | [pathtech.com.au](http://pathtech.com.au)

# LVL Starter Packs

Item	What's included	Small	Medium	Large
		(114-STP-S)	(114-STP-M)	(114-STP-L)
	1x Single Tube Reader (114-STR)	Yes	Yes	Yes
	1x SBS Format Rack Reader (114-DMTR-SBS)	No	Yes	Yes
	1x Single Tube Picker (114-CDC-1CH-XT-96) SBS 96	1	1	1
	1x Semi Automated Decapper (114-CDC-8CH-XT-4) 96	No	No	1
	Number of tube racks 300,500,1000ul XSX 96	10 racks	50 racks	100 racks

